



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Safe-Sense Coverall with Hood

**Product Code :**  
LBNY-00017-11000734



**Description :**

Sterile coveralls made of microporous polypropylene with polyethylene film feature elasticated ankles, hood, and back.

Sterile via gamma irradiation

PPE Category III

Zip front with sealable flap cover

Elasticated ankles, hood and back

Elastic cuffs with thumb loop

Disposable, full-body, one-piece protective clothing

Microporous polypropylene with polyethylene film

Low-linting and durable

Not made with natural rubber latex.

**Technical Specification :**

Safe-Sense Coverall with Hood



**Laboratory Instrument India**

