



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Light Grey Lab Coat with Exterior Pocket

**Product Code :**  
LBNY-00017-11000735



**Description :**

The left chest pocket with pencil stall and two lower packets will keep you stocked with pens and small tools while the one-piece, lined collar and button-front closure will keep you looking professional during rounds and time in the lab.

Primary Closure: Button

Pocket: Left chest exterior pocket, Bartacked pencil stall, Two lower exterior pockets

Non-yellowing buttons for long-lasting whiteness and professional appearance

Requires minimum ironing thanks to a wrinkle-resistant finish

Side vent openings for easy access to underneath attire

Fit Details: Fits true to size. Meant to be worn over clothes.

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Closure: Five buttons

Collar: One-piece, lined, notched lapel.

**Technical Specification :**

Light Grey Lab Coat with Exterior Pocket



**Laboratory Instrument India**

