



Email : sales@laboratoryinstrumentindia.com

Product Name :
Light Grey Lab Coat with Exterior Pocket

Product Code :
LBNY-00017-11000735



Description :

The left chest pocket with pencil stall and two lower packets will keep you stocked with pens and small tools while the one-piece, lined collar and button-front closure will keep you looking professional during rounds and time in the lab.

Primary Closure: Button

Pocket: Left chest exterior pocket, Bartacked pencil stall, Two lower exterior pockets

Non-yellowing buttons for long-lasting whiteness and professional appearance

Requires minimum ironing thanks to a wrinkle-resistant finish

Side vent openings for easy access to underneath attire

Fit Details: Fits true to size. Meant to be worn over clothes.

Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash

Closure: Five buttons

Collar: One-piece, lined, notched lapel.

Technical Specification :

Light Grey Lab Coat with Exterior Pocket



Laboratory Instrument India

