



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Onguard PVC Boot Covers

**Product Code :**  
LBNY-00017-11000736



**Description :**

Waterproof boot covers made of 0.2 mm thick PVC plastic uppers heat welded to 0.5 mm thick embossed PVC soles.

Waterproof

Color: Yellow.

15 in. high (38cm) boots with elastic top.

**Technical Specification :**

Onguard PVC Boot Covers



**Laboratory Instrument India**