



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens Lightweight Navy Long Sleeve Flame Resistant T-shirt

**Product Code :**  
LBNY-00017-11000737



**Description :**

Mens Lightweight Navy Long Sleeve Flame Resistant T-shirt is a CAT 2 AR/FR, compliant tee featuring a heat seal Bulwark triangle centered at the neck to show the wearer has on an FR undergarment.

Extra comfortable fabric

Underarm gusset for improved breathability

Rib-knit crewneck collar

Single spade chest pocket

Hazard: Flash Fire, Arc Flash

Wash Care: Home Wash

Fabric: Flame resistant, 7 oz. (235 g/m2) Interlock

Blend: 100% Cotton.

**Technical Specification :**

Mens Lightweight Navy Long Sleeve Flame Resistant T-shirt



**Laboratory Instrument India**

