



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Super Track Boot Covers

**Product Code :**  
LBNY-00017-11000739



**Description :**

For Cleanroom use  
With non-woven conductive/dissipative strip  
Elasticized seams and opening.  
Disposable  
Virgin polyethylene laminated to polypropylene  
High top  
Height: 7 in.

**Technical Specification :**

Super Track Boot Covers



**Laboratory Instrument India**