



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Ball and Ring

**Product Code :**  
LBNY-0003-1010033



**Description :**

Ball and Ring

**Technical Specification :**

It is used for demonstrating thermal expansion, comprising a captive brass ball secured to a mounted brass ring by a chain.

Ring mounted on rod with a wooden handle.

The ball passes through the ring when cold but will not pass through after being heated.

Brass ring with hard wood handle with a ball 25mm diameter attached by chain.



**Laboratory Instrument India**