



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Premium Coated Barrier Cleanroom Gowns

**Product Code :**  
LBNY-00017-11000745



**Description :**

Static dissipative, lint-free film coated gown for use in critical and low particulate ISO cleanroom class 6 environments.

Fabric is suitable for exposure to bloodborne pathogens, some hazardous liquids, and where liquid splash is a concern

Elastic wrists, elastic thumb loop, wrap-around neck ties, wrap-around waist ties

Latex free, lint-free film coated gown for use in critical and low particulate clean environments

Static dissipative

Fabric is constructed from a blend of polyethylene and polypropylene into a unique breathable membrane

Serged seam construction

Fabric meets or exceeds NIOSH, OSHA, ASHP and ONS guidelines for exposure to chemotherapy drugs.

**Technical Specification :**

Premium Coated Barrier Cleanroom Gowns



**Laboratory Instrument India**

