



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Class E Hi-Vis Insulated Bibs

**Product Code :**  
LBNY-00017-11000747



**Description :**

Made with tough 300D polyester oxford to withstand wear and tear  
PU transparent coating provides wind and water-resistant protection  
Scuff-resistant knee patches for extra durability  
Hi-Vis Insulated Bibs provide conspicuity, warmth, and durability for the toughest conditions.  
Adjustable straps and buckle release for comfort and fit  
1 zippered chest pocket and 2 side hand pockets for storage  
Boot zippers for easy on and off without removing footwear.  
ANSI Class E Complaint for low light conditions  
Black accents in high-wear areas to hide dirt and grime.

**Technical Specification :**

Class E Hi-Vis Insulated Bibs



**Laboratory Instrument India**