



Email : sales@laboratoryinstrumentindia.com

Product Name :
Class E Hi-Vis Insulated Bibs

Product Code :
LBNY-00017-11000747



Description :

Made with tough 300D polyester oxford to withstand wear and tear
PU transparent coating provides wind and water-resistant protection
Scuff-resistant knee patches for extra durability
Hi-Vis Insulated Bibs provide conspicuity, warmth, and durability for the toughest conditions.
Adjustable straps and buckle release for comfort and fit
1 zippered chest pocket and 2 side hand pockets for storage
Boot zippers for easy on and off without removing footwear.
ANSI Class E Complaint for low light conditions
Black accents in high-wear areas to hide dirt and grime.

Technical Specification :

Class E Hi-Vis Insulated Bibs



Laboratory Instrument India