



Email : sales@laboratoryinstrumentindia.com

Product Name :
Class E Thermal Pants

Product Code :
LBNY-00017-11000748



Description :

Class E Thermal Pants feature breathable Thinsulate warmth and waterproof comfort.
Elastic waistband with drawstring
Zippered leg openings for easy on/off
One right leg cargo pocket.
Weatherproof 100% polyester outer shell
Breathable PU coating
80 g Thinsulate lining
2 in. ANSI-compliant Reflective.

Technical Specification :

Class E Thermal Pants



Laboratory Instrument India