



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type R Class 2 Standard Vest

**Product Code :**  
LBNY-00017-11000753



**Description :**

Standard vest with choice of solid or mesh fabric and either hook and loop or zipper closure.

Material: Polyester.

Closure Type: Hook and Loop

Color: Lime

Size: 4X-Large to 5X-Large

**Technical Specification :**

Type R Class 2 Standard Vest



**Laboratory Instrument India**