



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Type R Class 2 Hi-Gloss Surveyors Vest

**Product Code :**  
LBNY-00017-11000754



**Description :**

Hi-Gloss, solid front, mesh back, zipper closure.

**Technical Specification :**

Front/mesh back  
Material: Polyester  
No. of Pockets: 2 inside, 4 outside.  
Color: Orange  
Size: 4X-Large to 5X-Large



**Laboratory Instrument India**