



Email : sales@laboratoryinstrumentindia.com

Product Name :
Traffic Officer Vest

Product Code :
LBNY-00017-11000755



Description :

Vest is fashioned from lime mesh with adjustable hook and loop sides
Adjustable hook and loop sides
Right chest mic tab
Polyester Lime Mesh
Zipper front closure
TRAFFIC OFFICER printed in reflective lettering on upper back.
2 in. silver Reflexite reflective material.

Technical Specification :

Traffic Officer Vest



Laboratory Instrument India