



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Traffic Officer Vest

**Product Code :**  
LBNY-00017-11000755



**Description :**

Vest is fashioned from lime mesh with adjustable hook and loop sides  
Adjustable hook and loop sides  
Right chest mic tab  
Polyester Lime Mesh  
Zipper front closure  
TRAFFIC OFFICER printed in reflective lettering on upper back.  
2 in. silver Reflexite reflective material.

**Technical Specification :**

Traffic Officer Vest



**Laboratory Instrument India**