



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Enhanced Visibility Multi Pocket Mesh Vest

**Product Code :**  
LBNY-00017-11000757



**Description :**

Created for maximum visibility.  
Zipper front closure reinforced with durable webbing  
Enhanced Visibility Multi Pocket Mesh Vest is made of 100% polyester mesh.  
Left and right chest mic tabs  
2 in. wide high performance reflective material with 3 in. high visibility contrasting lime  
Two lapel grommets.

**Pockets:**

Right chest 2-tier pencil pockets  
Left chest radio pocket (Style: Black/Lime); Left chest radio and clear Two lower outside cargo pockets with adjustable flaps (Style: B100 Black/Lime); Two lower outside cargo pockets with adjustable flaps and grommets for flagging tape (Style: B102 Royal Blue/Lime and B103 Red/Lime)  
ID pocket with flap (Style: Royal Blue/Lime and B103 Red/Lime)  
Two inside patch pocket.

**Technical Specification :**

Enhanced Visibility Multi Pocket Mesh Vest

---

# Laboratory Instrument India

