



Email : sales@laboratoryinstrumentindia.com

Product Name :
Enhanced Visibility Multi Pocket Mesh Vest

Product Code :
LBNY-00017-11000757



Description :

Created for maximum visibility.
Zipper front closure reinforced with durable webbing
Enhanced Visibility Multi Pocket Mesh Vest is made of 100% polyester mesh.
Left and right chest mic tabs
2 in. wide high performance reflective material with 3 in. high visibility contrasting lime
Two lapel grommets.

Pockets:

Right chest 2-tier pencil pockets
Left chest radio pocket (Style: Black/Lime); Left chest radio and clear Two lower outside cargo pockets with adjustable flaps (Style: B100 Black/Lime); Two lower outside cargo pockets with adjustable flaps and grommets for flagging tape (Style: B102 Royal Blue/Lime and B103 Red/Lime)
ID pocket with flap (Style: Royal Blue/Lime and B103 Red/Lime)
Two inside patch pocket.

Technical Specification :

Enhanced Visibility Multi Pocket Mesh Vest

Laboratory Instrument India

