



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Black Series Surveyors Vest

Product Code :
LBNY-00017-11000758



Description :

Premium Black Series Surveyors Vest with Class 2 Safety Compliance
Solid polyester front with mesh back
2 in. wide high performance reflective material laminated to 3 in. contrasting color
Heavy duty black bottom to help maintain a clean appearance
Zipper front closure reinforced with durable webbing
Left and right mic tabs (Size small only has left mic tab)
2 lapel grommets
Padded neck for comfort.

Technical Specification :

Premium Black Series Surveyors Vest



Laboratory Instrument India