



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Soft Scrubs SMS Long Sleeve Shirts

**Product Code :**  
LBNY-00017-11000759



**Description :**

Disposable long sleeve SMS shirts are comfortable to wear for long durations and are ideal for controlled environments.

Anti-static fabric is made of an opaque material that doesn't show through and resists abrasions

Convenient, fabric loop for easily hanging up clothing

Elastic wrists stretch to provide snug, secure fit

Round hemmed neck comfortably lays flat

Left chest pocket allows storage for eyeglasses, small instruments, and pens

Applications include pharmaceuticals, labs, healthcare, controlled environments, and more

Long sleeve SMS shirts offer comfort and safety to work long hours

Offer modesty with complete no show through coverage

Can easily be worn as a protective undergarment beneath a sterile gown

After use, quickly and easily dispose to reduce spread of germs, bioburdens, and other contagions.

**Technical Specification :**

Soft Scrubs SMS Long Sleeve Shirts



**Laboratory Instrument India**

