



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polypropylene Lab Coats, Three Pocket

Product Code :
LBNY-00017-11000760



Description :

Polypropylene lab coats protect from dust, dirt, grime, and other dry particulates in non-hazardous environments. Polypropylene lab coats feature three convenient pockets and provide economical, lightweight protection against dirt, grime, and particulates in non-hazardous environments.

Ideal for laboratories, pharmaceuticals, food processing, environmental cleanup sites, controlled environments

Breathable fabric helps prevent heat stress and keeps wearers cool over long hours

Three spacious pockets to carry eyeglasses, pens, and other small instruments

Low-cost protection.

Technical Specification :

Polypropylene Lab Coats, Three Pocket



Laboratory Instrument India