



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cleanroom Frocks

Product Code :
LBNY-00017-11000764



Description :

Ideal for non-hazardous work environments where reliable liquid and particulate protection is needed
Designed with tunnelized elastic to contain particle shed
Cleanroom frocks utilize state-of-the-art lamination technology to combine a soft non-woven fabric with a microporous film.
Front snap closures are easy to secure
Constructed with bound seams to eliminate raw edges and strengthen the seam
No pockets
Microporous, clean processed frocks available in individual or bulk packaging
Fabric keeps dust and liquids out while allowing air to flow freely
Bound mandarin collar for enhanced protection.

Technical Specification :

Cleanroom Frocks



Laboratory Instrument India

