



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cleanroom Frocks

**Product Code :**  
LBNY-00017-11000764



**Description :**

Ideal for non-hazardous work environments where reliable liquid and particulate protection is needed  
Designed with tunnelized elastic to contain particle shed  
Cleanroom frocks utilize state-of-the-art lamination technology to combine a soft non-woven fabric with a microporous film.  
Front snap closures are easy to secure  
Constructed with bound seams to eliminate raw edges and strengthen the seam  
No pockets  
Microporous, clean processed frocks available in individual or bulk packaging  
Fabric keeps dust and liquids out while allowing air to flow freely  
Bound mandarin collar for enhanced protection.

**Technical Specification :**

Cleanroom Frocks



**Laboratory Instrument India**

