



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coverall with Attached Hood

Product Code :
LBNY-00017-11000769



Description :

Provides Level C protective protection for a variety of chemical applications
Features inner sleeve and ankle cuffs, heavy-duty zipper, clear I.D. pocket, and reinforced crotch, elbows and knees
Constructed of 0.40mm nylon/polyester lining coated with heavy-duty PVC
Attached hood includes heavy-duty drawstring with cord lock.

Technical Specification :

Coverall with Attached Hood



Laboratory Instrument India