



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pleated Bouffant Cap

Product Code :
LBNY-00017-11000770



Description :

Lightweight and cool, breathable fabric
Elastic trimmed for a comfortable, snug fit
Disposable, full-coverage cleanroom bouffant cap.
Spunbonded polypropylene non-linting fabric
Elastic band is latex-free
Color: Sky blue
Opened diameter: 24 in. (60.96cm).

Technical Specification :

Pleated Bouffant Cap



Laboratory Instrument India