



Email : sales@laboratoryinstrumentindia.com

Product Name :
Polypropylene Fabric Grounding Ribbon Shoe Cover, 16 in.

Product Code :
LBNY-00017-11000772



Description :

Manufactured with 100% non-linting, non-woven spun bond polypropylene fabric. This shoe cover has stitched seams, sprayed latex anti-skid treads on the bottom for safety
It is commonly used for protecting against dirt, grime and certain dry particulates in nonhazardous environments and for protecting sensitive devices in the work area
All shoe covers are tested and manufactured in ISO compliant facilities under individual inspection and strict process control to ensure quality standard and product specifications
This product is recommended for use in a cleanroom Class 10,000+ (ISO 7+) critical environment
19.6 in. long conductive ribbon stitched on a solid seam to provide ESD grounding
Lightweight, strong, and durable.

Technical Specification :

Polypropylene Fabric Grounding Ribbon Shoe Cover, 16 in.



Laboratory Instrument India