



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Heavyweight Polyethylene Shoe Covers

**Product Code :**  
LBNY-00017-11000773



**Description :**

Disposable cross-linked polyethylene shoe covers are manufactured with 100% cross-linked polyethylene with heat-sealed seams

Shoe covers are tested and manufactured in ISO compliant facilities under inspection and stringent process control to ensure quality standard and product specification compliance

It is the most economical shoe cover available that provides full liquid protection and is commonly used in microelectronics and pharmaceutical industries

Heavyweight disposable shoe covers are skillfully crafted from 100% cross-linked polyethylene and incorporate heat-sealed seams. Strong, durable, impervious to liquids, and suitable for ISO 5-6 (Class 100-1,000) critical environments.

This product is recommended for use in a cleanroom Class 100 -1.000 (ISO 5-6) critical environment.

**Technical Specification :**

Heavyweight Polyethylene Shoe Covers



**Laboratory Instrument India**

