



Email : sales@laboratoryinstrumentindia.com

Product Name :
Front-Load Dispensers with Opening for Access

Product Code :
LBNY-00017-11000774



Description :

Material: Acrylic
Product Style: Front Loading with Opening
Color: Clear
Depth (Metric): 55.9 cm
Depth (English): 22 in.
Width (Metric): 30.5 cm
Width (English): 10.5 in.
Height (Metric): 40.6 cm
For Use With (Application): Designed to fit Metro rack.

Technical Specification :

Front-Load Dispensers with Opening for Access



Laboratory Instrument India