




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Tians Epic Shoe Cover

**Product Code :**  
LBNY-00017-11000776



**Description :**

Made from a durable soild polyethylene film  
Fabric is constructed from a proprietary polyethylene blend, and is embossed for additional skid resistance  
Clean and economical non-linting shoe cover for critical environments  
Speeds up weighing process  
Fabric is suitable for exposure to bloodborne pathogens, some hazardous liquids, and where liquid splash is a concern  
Welded seams  
Lint free disposable cleanroom shoe cover with heat sealed seams  
Seamless bottom  
For use in critical and low particulate clean environments.

**Technical Specification :**

Tians Epic Shoe Cover



**Laboratory Instrument India**

