




Email : sales@laboratoryinstrumentindia.com

Product Name :
Tians Epic Shoe Cover

Product Code :
LBNY-00017-11000776



Description :

Made from a durable soild polyethylene film
Fabric is constructed from a proprietary polyethylene blend, and is embossed for additional skid resistance
Clean and economical non-linting shoe cover for critical environments
Speeds up weighing process
Fabric is suitable for exposure to bloodborne pathogens, some hazardous liquids, and where liquid splash is a concern
Welded seams
Lint free disposable cleanroom shoe cover with heat sealed seams
Seamless bottom
For use in critical and low particulate clean environments.

Technical Specification :

Tians Epic Shoe Cover



Laboratory Instrument India

