



Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium Coated Barrier Cleanroom Gowns

Product Code :
LBNY-00017-11000777



Description :

Fabric is suitable for exposure to bloodborne pathogens, some hazardous liquids, and where liquid splash is a concern

Elastic wrists, elastic thumb loop, wrap-around neck ties, wrap-around waist ties

Latex free, lint-free film coated gown for use in critical and low particulate clean environments

Static dissipative, lint-free film coated gown for use in critical and low particulate ISO cleanroom class 6 environments.

Fabric is constructed from a blend of polyethylene and polypropylene into a unique breathable membrane

Static dissipative

Serged seam construction

Fabric meets or exceeds NIOSH, OSHA, ASHP and ONS guidelines for exposure to chemotherapy drugs.

Technical Specification :

Premium Coated Barrier Cleanroom Gowns



Laboratory Instrument India

