



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
65 Cal Bib Overalls

**Product Code :**  
LBNY-00017-11000778



**Description :**

Fiber Provides lightweight, durable, heat and flame-resistant solutions. Made of 34% aramid, 33% lyocell, 31% modacrylic, and 2% anti-static materials is engineered to absorb moisture and dry quickly

Fabric/Material: Multilayer

FR fabric will not melt, drip, or ignite

Quick disconnect and easy adjustable shoulder straps for quick entry

Expandable leg opening to fit over boots

FR hook and loop calf closure for easy boot entry

Quick disconnect buckles at each shoulder

Arc Rating: 65 cal/cm<sup>2</sup>

CAT 4

Inherently flame resistant for FR protection

Garments embroidered with arc rating

Available in a convenient kit.

**Technical Specification :**

65 Cal Bib Overalls

---

# Laboratory Instrument India

