



Email : sales@laboratoryinstrumentindia.com

Product Name :
40 Cal Bib Overalls SDP

Product Code :
LBNY-00017-11000780



Description :

FR fabric will not melt, drip, or ignite
Quick disconnect and easy adjustable shoulder straps
Heavy duty 26 in. zipper closure on each leg
Garments embroidered with arc rating
Available in a convenient kit
National Safety Apparel are 50% lighter than standard 40 cal/cm² Arc Flash PPE
Inherently flame resistant for FR protection
Expandable leg opening to fit over boots

Technical Specification :

40 Cal Bib Overalls SDP



Laboratory Instrument India