



Email : sales@laboratoryinstrumentindia.com

Product Name :
FR Fleece 3-in-1 Hood

Product Code :
LBNY-00017-11000783



Description :

Designed for cold weather conditions, the FR Fleece 3-in-1 Hood combines a fleece hood with an adjustable inner facial liner and neck gaiter
Fabric/Material: 8.5 oz. IIIA Fleece (Hood) / 5.5 oz. FR Control 2.0 (Neck gaiter/face cover)
Adjustable drawstring hood
Combines a fleece hood with an adjustable inner facial liner/neck gaiter
Breathable and quick drying FR Fabric
Wicks moisture away from the skin for comfort and warmth (inner liner)
Designed for cold weather conditions
Several unique ways to wear
Arc Rating: 24 cal/cm²
CAT 2
Warmth without the weight.

Technical Specification :

FR Fleece 3-in-1 Hood



Laboratory Instrument India

