



Email : sales@laboratoryinstrumentindia.com

Product Name :
Performance Arc Flash Kit, 40 Cal

Product Code :
LBNY-00017-11000784



Description :

Using comfort science to keep the wearer cool, dry, comfortable, and unrestricted while working
Ergonomic design for ease of movement and flexibility
Single-layer 40 cal arc flash protection
Up to 20% lighter, 48% thinner, and 56% more flexible than other premium arc flash suits
26 in. leg zippers for easy donning & doffing
Enhanced visibility with FR silver segmented reflective trim
Face shield technology features a greater range of visibility
Gloves are not included.
Performance Arc Flash Kit, 40 Cal, CAT 4, Includes Jacket, Bib, Lift Front Hood, Backpack.

Technical Specification :

Performance Arc Flash Kit, 40 Cal



Laboratory Instrument India