



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Navy Lab Coat with Exterior Pocket

Product Code :
LBNY-00017-11000785



Description :

Mens Navy Lab Coat with Exterior Pocket and Four-Button Closure
Blend: 80% Polyester / 20% Combed Cotton
Finish: Wrinkle Resistant
Closure: Five buttons
Collar: One-piece, lined, notched lapel
Primary Closure: Button
Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash
Fabric: 5 oz. Performance Blend Poplin
Length: RG - 41.5 in.; LN - 45.5 in. (all sizes - back)
One chest pocket and two lower pockets provide storage for tools and personal items
Side vent openings for easy access to underneath attire
Pocket: Left chest exterior pocket, Bartacked pencil stall, Two lower exterior pockets
Non-yellowing buttons for long-lasting whiteness and professional appearance
Requires minimum ironing thanks to a wrinkle-resistant finish
41.5 in. back length on all regular sizes (45.5 in. on longs).

Technical Specification :

Mens Navy Lab Coat with Exterior Pocket

Laboratory Instrument India

