



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Mens Navy Lab Coat with Exterior Pocket

**Product Code :**  
LBNY-00017-11000785



**Description :**

Mens Navy Lab Coat with Exterior Pocket and Four-Button Closure  
Blend: 80% Polyester / 20% Combed Cotton  
Finish: Wrinkle Resistant  
Closure: Five buttons  
Collar: One-piece, lined, notched lapel  
Primary Closure: Button  
Wash Care: Industrial Laundry - Heavy Soil, Industrial Laundry - Light Soil, Home Wash  
Fabric: 5 oz. Performance Blend Poplin  
Length: RG - 41.5 in.; LN - 45.5 in. (all sizes - back)  
One chest pocket and two lower pockets provide storage for tools and personal items  
Side vent openings for easy access to underneath attire  
Pocket: Left chest exterior pocket, Bartacked pencil stall, Two lower exterior pockets  
Non-yellowing buttons for long-lasting whiteness and professional appearance  
Requires minimum ironing thanks to a wrinkle-resistant finish  
41.5 in. back length on all regular sizes (45.5 in. on longs).

**Technical Specification :**

Mens Navy Lab Coat with Exterior Pocket

---

# Laboratory Instrument India

