



Email : sales@laboratoryinstrumentindia.com

Product Name :
Revolite Coat 40 Cal 49 in. Inseam Long

Product Code :
LBNY-00017-11000786



Description :

40 cal flash coat that is 50% lighter than our standard 40 cal/cm² arc flash kits
FR hook & loop front closure with drawstring waist
Stand up collar for extra coverage, that can be turned down
Knit cuffs for a comfortable, secure fit
Garments embroidered with arc rating
FR fabric will not melt, drip, or ignite
Arc Rating: 40 cal/cm²
CAT 4
Fabric/Material: 11.2 oz.
50% lighter than our standard 40 cal/cm² arc flash kits
Available in a convenient kit

Technical Specification :

Revolite Coat 40 Cal 49 in. Inseam Long



Laboratory Instrument India

