



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Revolite Coat 40 Cal 49 in. Inseam Long

**Product Code :**  
LBNY-00017-11000786



**Description :**

40 cal flash coat that is 50% lighter than our standard 40 cal/cm<sup>2</sup> arc flash kits  
FR hook & loop front closure with drawstring waist  
Stand up collar for extra coverage, that can be turned down  
Knit cuffs for a comfortable, secure fit  
Garments embroidered with arc rating  
FR fabric will not melt, drip, or ignite  
Arc Rating: 40 cal/cm<sup>2</sup>  
CAT 4  
Fabric/Material: 11.2 oz.  
50% lighter than our standard 40 cal/cm<sup>2</sup> arc flash kits  
Available in a convenient kit

**Technical Specification :**

Revolite Coat 40 Cal 49 in. Inseam Long



**Laboratory Instrument India**

