



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Contractor 40 Cal Hood With Hat 9 oz. Khaki 2 Ply 88/12

**Product Code :**

LBNY-00017-11000791



**Description :**

Hood with slotted adapted and breathable FR fabric  
Anti-fog faceshield  
Slotted adapter  
Breathable FR fabric  
Inherently flame resistant  
One size fits all  
Available in a convenient kit.

**Technical Specification :**

Contractor 40 Cal Hood With Hat 9 oz. Khaki 2 Ply 88/12



**Laboratory Instrument India**