



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Dark Brown Industrial Work Pant

Product Code :
LBNY-00017-11000793



Description :

The comprised of heavyweight 8.5 oz. twill.
Blend: Body: 65% Polyester / 35% Cotton.
Work pants built of sturdy fabric for hard wear and easy care. Tradesmen swear by the comfort and durability.
Finish: Wrinkle Resistant, Stain Release
Closure: Metal hook and eye
Pocket: Welt back pockets
Number of Pockets: 4
Primary Closure: Hook and Eye
Wash Care: Home Wash
Distinctive tunnel belt loops
Product Collection: Flex
Fabric: Body: 8.5 oz. Twill
Waistband: Inner-lined, 2 in.
Gender: Male.

Technical Specification :

Mens Dark Brown Industrial Work Pant

Laboratory Instrument India

