



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Mens White Long Sleeve Industrial Work Shirt

**Product Code :**

LBNY-00017-11000794



**Description :**

Easy-care, permanent-press polyester/cotton poplin resists wrinkling and fading. The ideal uniform shirt made to withstand the rigors of any job and industrial laundering.

Work shirt that is rugged, tough, comfortable, and presents a sharp appearance.

Wash Care: Industrial Laundry - Heavy Soil

Fabric: Body: 4.25 oz. Poplin

Blend: Body: 65% Polyester / 35% Cotton

Finish: Wrinkle Resistant

Number of Pockets: 2

Primary Closure: Button

Collar: Two-piece convertible collar with sewn in stays

Gender: Male

Pocket: Hex chest pockets with pencil stall

Country of Origin: Imported, Imported Fabric.

**Technical Specification :**

Mens White Long Sleeve Industrial Work Shirt

---

# Laboratory Instrument India

