



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Desert Sand Industrial Work Pant

Product Code :
LBNY-00017-11000795



Description :

Work pants built of sturdy fabric for hard wear and easy care. Tradesmen swear by the comfort and durability. The comprised of heavyweight 8.5 oz. twill.
Fabric: Body: 8.5 oz. Twill
Blend: Body: 65% Polyester / 35% Cotton
Finish: Wrinkle Resistant, Stain Release
Closure: Metal hook and eye
Pocket: Welt back pockets
Number of Pockets: 4
Distinctive tunnel belt loops
Product Collection: Flex
Primary Closure: Hook and Eye
Wash Care: Home Wash
Waistband: Inner-lined, 2 in.
Country of Origin: Imported, Imported Fabric
Gender: Male.

Technical Specification :

Mens Desert Sand Industrial Work Pant

Laboratory Instrument India

