



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens Dark Navy Premium Industrial Multi-Use Pocket Pant

Product Code :
LBNY-00017-11000798



Description :

Multi-use pocket pants offer additional pocket storage with a more professional look. Designed for durability, on the job and in industrial laundries.

Fabric: Body: 7.25 oz. Twill

Blend: Body: 65% Polyester / 35% Cotton

Finish: Wrinkle Resistant, Moisture Wicking, Stain Release

Closure: Metal tack button, Heavy duty brass ratcheting sipper won't break under stress

Number of Pockets: 5

Primary Closure: Button

Wash Care: Industrial Laundry - Heavy Soil

Pocket: Reinforced front pocket bags, Welt back pockets, Multi-use side pocket

Waistband: Inner-lined, 1.5 in.

E-Z Touch

It's a fade-resistant, permanent-press, polyester/cotton twill with Stay-Dark Technology.

Country of Origin: Imported, Imported Fabric

Fit Details: Sits slightly below waist, Extra room in seat and thigh, Straight leg

Gender: Male.

Technical Specification :

Mens Dark Navy Premium Industrial Multi-Use Pocket Pant



Laboratory Instrument India