



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Womens Navy Industrial Work Pant

**Product Code :**  
LBNY-00017-11000799



**Description :**

The tough polyester/cotton twill resists wrinkles and releases stains. It is a clean, straightforward look versatile enough for women in all types of work.

Work pants that sit slightly below the waist with a traditional fit through the seat and thigh.

Fabric: Body: 8.5 oz. Twill

Blend: Body: 65% Polyester / 35% Cotton

Finish: Wrinkle Resistant, Stain Release

Closure: Metal hook and eye, Ratcheting brass zipper

Number of Pockets: 4

Primary Closure: Hook and Eye

Wash Care: Home Wash

Resists wrinkles

Stain release

Pocket: Welt back pockets, Dyed to match button secures left back pocket

Waistband: Inner-lined, two-piece waistband, 2 in.

Two back welt pockets - left with button closure

Distinctive tunnel belt loops

Country of Origin: Imported, Imported Fabric

Fit Details: Sits slightly below the waist, Traditional fit in seat and thigh

Gender: Female.

**Technical Specification :**

Womens Navy Industrial Work Pant



**Laboratory Instrument India**