



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mens Silver Industrial Work Pant

**Product Code :**  
LBNY-00017-11000800



**Description :**

Work pants built of sturdy fabric for hard wear and easy care. Tradesmen swear by the 874's comfort and durability.

Fabric: Body: 8.5 oz. Twill

Blend: Body: 65% Polyester / 35% Cotton

Finish: Wrinkle Resistant, Stain Release

Closure: Metal hook and eye

The comprised of heavyweight 8.5 oz. twill.

Waistband: Inner-lined, 2 in.

Distinctive tunnel belt loops

Product Collection: Flex

Number of Pockets: 4

Primary Closure: Hook and Eye

Wash Care: Home Wash

Pocket: Welt back pockets

Country of Origin: Imported, Imported Fabric

Gender: Male.

**Technical Specification :**

Mens Silver Industrial Work Pant



**Laboratory Instrument India**