




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Premium FR Phase Change Cooling Vest with Packs

**Product Code :**  
LBNY-00017-11000801



**Description :**

Premium Phase Change Cooling Vest with Packs keeps workers cool and safe while maintaining a cool temperature for up to 4 hours when activated.  
Features two elastic pen loops on chest.  
Vest made with modacrylic cotton blend FR fabric  
Flexible fit  
Adjustable hook and loop closures on both shoulders and waist  
Extended vest surface area enhances cooling power  
Quick recharge in a little as 10 minutes  
Honeycomb design eliminates stiffness  
Flame resistant  
Immediate and long-lasting cooling relief  
Packs maintain constant 64°F (18°C) for up to 4 hours  
Designed for more visibility on busy streets or job sites.

**Technical Specification :**

Premium FR Phase Change Cooling Vest with Packs

---

# Laboratory Instrument India

