




Email : sales@laboratoryinstrumentindia.com

Product Name :
Premium FR Phase Change Cooling Vest with Packs

Product Code :
LBNY-00017-11000801



Description :

Premium Phase Change Cooling Vest with Packs keeps workers cool and safe while maintaining a cool temperature for up to 4 hours when activated.
Features two elastic pen loops on chest.
Vest made with modacrylic cotton blend FR fabric
Flexible fit
Adjustable hook and loop closures on both shoulders and waist
Extended vest surface area enhances cooling power
Quick recharge in a little as 10 minutes
Honeycomb design eliminates stiffness
Flame resistant
Immediate and long-lasting cooling relief
Packs maintain constant 64°F (18°C) for up to 4 hours
Designed for more visibility on busy streets or job sites.

Technical Specification :

Premium FR Phase Change Cooling Vest with Packs

Laboratory Instrument India

