



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Anti-Skid Polypropylene Shoe Covers

**Product Code :**  
LBNY-00017-11000803



**Description :**

Manufactured with 100% non-linting, non-woven spun bond polypropylene fabric with stitched seams and sprayed latex anti-skid treads on the sole.  
Recommended for use in a cleanroom Class 10,000+ (ISO 7+) critical environment  
Commonly used for protecting against dirt, grime, and certain dry particulates in nonhazardous environments.  
All shoe covers are tested and manufactured in ISO compliant facilities under individual inspection and strict process control to ensure quality standard and product specifications.

**Technical Specification :**

Anti-Skid Polypropylene Shoe Covers



**Laboratory Instrument India**