



Email : sales@laboratoryinstrumentindia.com

Product Name :
Microporous Lab Coats with Full Zipper

Product Code :
LBNY-00017-11000804



Description :

Lab coats are manufactured with non-linting microporous film bonded with a non-woven fabric to provide barrier protection against non-toxic particulates, liquids, sprays, and dust.

Ideal for use in a cleanroom Class 1-10 critical environment that requires high contamination control and cleanliness.

Tested and manufactured in ISO compliant facilities to ensure quality standard and product specification compliance

Recommended for either stringent industrial or demanding cleanroom applications.

Highly breathable and comfortable

Knee length with zipper closure, one pocket, elastic cuff, and a single-ply shirt collar design.

Technical Specification :

Microporous Lab Coats with Full Zipper



Laboratory Instrument India