



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Weight Polyethylene Shoe Cover

Product Code :
LBNY-00017-11000805



Description :

Provides full liquid protection and is commonly used in microelectronics and pharmaceutical industries. Manufactured with 100% cross-linked polyethylene with heat-sealed seams. Strong, durable, impervious to non-hazardous liquids, and critical environment compatible.

Technical Specification :

Standard Weight Polyethylene Shoe Cover



Laboratory Instrument India